

Fabrication of Single-walled Carbon Nanotube (SWNT) Sorbent for Volatile Organic Compounds (VOCs) Sampling

Jonghwa Oh^{*1}, Claudiu T. Lungu¹, Evan Floyd²

* Presenter

¹ University of Alabama at Birmingham

² University of Oklahoma

INTRODUCTION: VOCs Sampling and Analysis



SKC Anasorb CSC Sorbent Tube

- Sampling (Adsorption)
 - (Ad)sorbent: Charcoal or Activated Carbon
 - High adsorption capacity & affordable
- Analysis (Desorption and Measurement)
 - Chemical extraction
 - Long desorption time (i.e., > 30 min) in a toxic solvent
 - Somewhat expensive GC (gas chromatography) vials
 - Low sensitivity
 - Thermal extraction
 - Desorption directly in GC system; much more sensitive
 - Reliability issue (e.g., poor sealing)
 - Expensive system

INTRODUCTION: New Desorption Technique



- PhotoThermal Desorption (PTD)
 - Visible light (photo flash) thermally desorbing an analyte & quantified by photoionization detector
 - An alternative for VOCs analysis of *passive sampler*
 - Elimination of sample preparation time and cost
 - Bridging the sensitivity gap between chemical desorption and thermal desorption



10th IOHA International
Scientific Conference

IOHA London 2015, 27-30 April
2015, Hilton London Metropole




INTRODUCTION: Potential Sorbent



Carbon Nanotubes (CNTs)

- *Large surface area* (300-1,000 m²/g)
 - Tunable pore diameter by fabrication process
- *Superior thermal, electrical, optical conductivity*
 - A wide variety of applications (e.g., mechanical, structural, thermal, electrical, biomedical, etc.)
- *Single-walled Carbon Nanotubes (SWNTs)*
 - Considered as prime material for gas sorption
 - Higher in thermal/electrical conductivity (3000 W/m·K)¹



One, Two, and Three walls of CNTs

¹ Thermal conductivity of gold: 310 W/m·K

PURPOSE

- Develop a CNT sorbent coupled with photothermal desorption technique to be used in a VOC passive sampler
 - Investigate fabrication methods of SWNTs to a self-standing form (i.e., buckypaper)
 - Examine adsorption efficiencies (i.e., surface area, pore size, and toluene adsorption isotherm)

METHODS:

Fabrication of Buckypaper (BP)

- Arc Discharge (AD) SWNT solution¹
- Suspension and filtration
 - Suspension with solvents/surfactants (i.e., acetone and methanol)
 - Vacuum filtration
 - Additional cleaning with deionized water and solvents²
- Three types of BP samples (50 mg each, n=2)
 - Non-cleaned (acetone suspension³ without cleaning)
 - Acetone-cleaned (acetone suspension³ with cleaning)
 - Methanol-cleaned (methanol suspension³ with cleaning)

¹ Suspended in 1 % surfactants (i.e., sodium cholate and sodium dodecyl sulfate)


² Repeated 2 times with 250 mL deionized water and 50 mL solvent

³ 15 hrs of suspension in 400 mL solvent

METHODS

SWNTs
suspended
in solvent

Membrane
filter



Vacuum

METHODS:

Adsorption Efficiency


- Surface area & Pore size
 - ASAP2020 Physisorption Analyzer
 - Nitrogen as adsorbate
 - 77 K (-196 °C)
- Toluene Adsorption Isotherm¹
 - Physisorption Analyzer
 - Toluene as adsorbate
 - 25 °C



¹ Analyzed at Micromeritics® because of a special set-up required

RESULTS: Fabrication of Buckypaper

◦ Non-cleaned BP



◦ Acetone-cleaned BP



◦ Methanol-cleaned BP



RESULTS: Adsorption Efficiency



- Surface area¹ and Pore size of Arc Discharge² BP

| Fabrication Method | NOT-CLEANED BP (Acetone suspension & Filtration w/o Cleaning) | ACETONE-CLEANED BP (Acetone suspension & Filtration with Cleaning) | METHANOL- CLEANED BP (Methanol suspension & Filtration with Cleaning) |
|---|--|---|--|
| Surface Area (m ² /g) ³ | 45 ± 2 | 217 ± 27 (349 ± 10)⁴ | 348 ± 13 (421 ± 6)⁴ |
| Pore Size (nm) | 14.9 ± 0.0 | 9.1 ± 0.4 (9.8 ± 0.1)⁴ | 8.0 ± 0.1 (8.6 ± 0.1)⁴ |

- ¹ Surface area measurement was repeated 2 times per sample and averaged
- ² Suspended in 1 % surfactants (i.e., sodium cholate and sodium dodecyl sulfate)
- ³ Determined by Brunauer, Emmett and Teller (BET) theory
- ⁴ Measured after toluene adsorption isotherm in which 960 min of degassing was performed at 300 °C



10th IOHA International
Scientific Conference



IOHA London 2015, 27-30 April
2015, Hilton London Metropole



RESULTS: Adsorption Efficiency

- Toluene Adsorption Isotherm¹
 - Acetone-cleaned BP (150 cm³/g STP)

- Methanol-cleaned BP (187 cm³/g STP)



¹ Measurement was performed with one representative sample

RESULTS: Adsorption Efficiency



- Toluene Adsorption Capacity¹

| Fabrication Method | Not-Cleaned BP | Acetone-Cleaned BP | Methanol-Cleaned BP |
|--|----------------|--------------------|---------------------|
| Adsorption Quantity (cm ³ toluene /g BP at STP) | — ² | 150 | 187 |
| Adsorption Capacity ³ (cm ³ toluene / g BP) | — ² | 0.72 | 0.89 |
| Adsorption Capacity (mg toluene / g BP) | — ² | 618 | 770 |

- ¹ Measurement was repeated 2 times and averaged
- ² Data was not included because of considerably lower surface area
- ³ Conversion of gas volume at STP to liquid volume at analytical temperature using a toluene density conversion factor of 0.0047677



10th IOHA International
Scientific Conference

IOHA London 2015, 27-30 April
2015, Hilton London Metropole



Conclusion/ Future Work

- Cleaning process was necessary for this material (i.e., arc discharge SWNT solution) and methanol cleaning fabrication method was the most efficient
- Toluene adsorption capacity was proportional to surface area
- Increase surface area by fabrication (e.g., annealing)
- Examine other types of sorbents (e.g., HiPco (high pressure CO) SWNTs, carbonaceous particles, etc.)
- Test PhotoThermal Desorption (PTD) efficiency of the samples
- Build a passive sampler unit for photothermal desorption

LONDON 2015



감사합니다 Natick
Grazie Danke Ευχαριστίες Dalu
Thank You Köszönöm
Спасибо Dank Gracias
谢谢 Merci Seé
ありがとう

Obrigado

LONDON 2015



10th IOHA International
Scientific Conference

IOHA London 2015, 27-30 April
2015, Hilton London Metropole



The Chartered Society for
Worker Health Protection



10th IOHA International Scientific Conference

LONDON 2015



April 25th-30th 2015
London Metropole Hotel



IOHA & BOHS 2015
London: Building on Occupational Hygiene Together

25-30 April 2015
London Metropole Hotel
www.iohalondon2015.org